

Anti-Tle1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q04724

Item Number: CM0008832



Introduction

High-quality rabbit monoclonal TLE1 antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Suitable for gene expression, cell signaling, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TLE1 Rabbit mAb
Remarks/Alias	ESG; ESG1; GRG1; TLE-1; TLE1
Species	Human
GeneID (Human)	7088
GeneID	7088
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human TLE1 (Q04724)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04724
Protein Weight	83 kDa
Shipping	Ice bag