

Anti-Tle4 Monoclonal Antibody For Wb And Elisa - Q04727

Item Number: CM0007485

Introduction

High-quality rabbit monoclonal antibody against human TLE4, reactive with human and mouse. Validated for WB and ELISA, ideal for transcriptional repression and Wnt signaling research.

[Learn More](#)

Parameter	Specification
Product Name	TLE4 Rabbit mAb
Remarks/Alias	ESG, BCE1, ESG4, GRG4, BCE-1, Grg-4, E(spl), E(spl), TLE4
Species	Human
Gene ID	7091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human TLE4 (NP_001269677)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04727
Protein Weight	84kDa
Shipping	Ice bag