

Anti-Tle4 Rabbit Polyclonal Antibody For Wb, Elisa - Q04727

Item Number: CM0022869

Introduction

High-quality TLE4 rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to rat. Supplied as a rabbit polyclonal antibody, it is suitable for detecting human TLE4 protein in studies of transcriptional repression in Wnt signaling and development. - Q04727

[Learn More](#)

Parameter	Value
Product Name	TLE4 Rabbit pAb
Remarks/Alias	ESG, BCE1, ESG4, GRG4, BCE-1, Grg-4, E(spl), E(spl), TLE4
Species	Human
GeneID (Human)	7091
GeneID	7091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 173-243 of human TLE4 (NP_001269677)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q04727
Protein Weight	84kDa
Shipping	Ice bag