

Anti-Tlk1 Polyclonal Antibody For Wb, Elisa - Q9Uki8

Item Number: CM0024079

Introduction

Rabbit polyclonal antibody targeting TLK1 (PKU-beta), suitable for WB and ELISA. Cross-reacts with human, mouse, and rat TLK1. TLK1 is a kinase involved in chromatin assembly and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	TLK1 Rabbit pAb
Remarks/Alias	PKU-beta; TLK1
Species	Human
GeneID	9874
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human TLK1 (NP_036422.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKI8
Protein Weight	87kDa
Shipping	Ice bag