

Anti-Tlk1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9UKi8

Item Number: CM0013023

Introduction

Rabbit polyclonal antibody targeting human TLK1 (PKU-beta), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for cell cycle, DNA damage repair, and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	TLK1 Rabbit pAb
Remarks/Alias	PKU-beta; TLK1
Species	Human
Gene ID (Human)	9874
Gene ID	9874
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-230 of human TLK1 (NP_036422.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKI8
Protein Weight	87kDa
Shipping	Ice bag