

Anti-Tlr1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15399

Item Number: CM0015943



Introduction

Rabbit polyclonal antibody against human TLR1, suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for innate immunity and inflammatory response research.

[Learn More](#)

Parameter	Value
Product Name	TLR1 Rabbit pAb
Remarks/Alias	TIL; CD281; rsc786; TIL. LPRS5; TLR1
Species	Human
Gene ID (Human)	7096
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-290 of human TLR1 (NP_003254.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15399
Protein Weight	90kDa
Shipping	Ice bag