

Anti-Tlr2 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - O60603

Item Number: CM0012094



Introduction

High-quality anti-TLR2 rabbit monoclonal antibody validated for WB, IF/ICC, IF-P, and ELISA. Suitable for human and mouse samples. Ideal for innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	TLR2 Rabbit mAb
Remarks/Alias	TIL4; CD282
Species	Human
GeneID	7097
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 290-390 of human TLR2 (NP_003255.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60603
Protein Weight	90kDa
Shipping	Ice bag