

Anti-Tlr2 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - O60603

Item Number: CM0002928



Introduction

High-quality TLR2 Rabbit Monoclonal Antibody (O60603) for WB, IF/ICC, IF-P, and ELISA applications, cross-reactive with human and mouse. Suitable for studying innate immunity, TLR2 signaling, and bacterial infection pathways.

[Learn More](#)

Parameter	Value
Product Name	TLR2 Rabbit mAb
Remarks/Alias	TIL4; CD282; TLR2
Species	Human
Gene ID	7097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-784 of human TLR2 (NP_003255.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60603
Protein Weight	90kDa
Shipping	Ice bag