


Anti-Tlr3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15455

Item Number: CM0017639



Introduction

TLR3 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. It recognizes Human, Mouse, Rat TLR3 and is ideal for innate immunity and viral infection research.

[Learn More](#)

Parameter	Description
Product Name	TLR3 Rabbit pAb
Remarks/Alias	CD283; IIAE2; IMD83; TLR3
Species	Human
Gene ID (Human)	7098
Gene ID	7098
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human TLR3 (NP_003256.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15455
Protein Weight	104kDa
Shipping	Ice bag