

Anti-Tlr4 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - 000206

Item Number: CM0018147



Introduction

TLR4 Rabbit pAb is a rabbit polyclonal antibody against human TLR4, validated for WB, IF/ICC, IP, and ELISA. Reacts with human and mouse TLR4; ideal for studying innate immunity and inflammatory signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	TLR4 Rabbit pAb
Remarks/Alias	TOLL; CD284; TLR-4; ARMD10; TLR4
Species	Human
GeneID (Human)	7099
GeneID	7099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-228 of human TLR4 (NP_612564.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	000206
Protein Weight	96kDa
Shipping	Ice bag