


Anti-Tlr5 Rabbit Monoclonal Antibody For Wb And Elisa - O60602

Item Number: CM0008840



Introduction

High-quality rabbit monoclonal antibody against human TLR5, validated for Western blot (WB) and ELISA. Recognizes human TLR5 protein (~98 kDa). Suitable for innate immunity and inflammatory disease research. O60602.

[Learn More](#)

Parameter	Value
Product Name	TLR5 Rabbit mAb
Remarks/Alias	SLE1; TIL3; SLEB1; MELIOS; TLR5
Species	Human
Gene ID	7100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 258-404 of human TLR5 (O60602)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60602
Protein Weight	98 kDa
Shipping	Ice bag