


Anti-Tlr5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60602

Item Number: CM0021751



Introduction

Anti-TLR5 rabbit polyclonal antibody targeting the human toll-like receptor 5 (O60602). Validated for WB, IHC-P, and ELISA, this antibody aids in innate immunity and inflammatory response research.

[Learn More](#)

Parameter	Value
Product Name	TLR5 Rabbit pAb
Remarks/Alias	SLE1; TIL3; SLEB1; MELIOS; TLR5
Species	Human
GeneID (Human)	7100
GeneID	7100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 714-858 of human TLR5 (NP_003259.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O60602
Protein Weight	98kDa
Shipping	Ice bag