


Anti-Tlr6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2C9

Item Number: CM0025736



Introduction

TLR6 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Recognizes endogenous TLR6 in human and mouse samples. Ideal for innate immunity and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	TLR6 Rabbit pAb
Remarks/Alias	CD286; TLR6
Species	Human
Gene ID (Human)	10333
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 697-796 of human TLR6 (NP_006059.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2C9
Protein Weight	92kDa
Shipping	Ice bag