


Anti-Tlr7 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9Nyk1

Item Number: CM0013275



Introduction

Rabbit polyclonal antibody against human TLR7 (SWISS: Q9NYK1). Validated for WB, IF/ICC, IF-P, and ELISA. Reacts with human, mouse, and rat TLR7. Immunogen corresponds to the extracellular domain (aa 30-230). Ideal for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	TLR7 Rabbit pAb
Remarks/Alias	IMD74; SLEB17; TLR7-like; TLR7
Species	Human
GeneID (Human)	51284
GeneID	51284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-230 of human TLR7 (NP_057646.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYK1
Protein Weight	121kDa
Shipping	Ice bag