

Anti-Tlr9 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nr96

Item Number: CM0017129



Introduction

Rabbit polyclonal antibody targeting human Toll-like receptor 9 (TLR9/CD289). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for innate immunity research.

[Learn More](#)

Parameter	Specification
Product Name	TLR9 Rabbit pAb
Remarks/Alias	CD289; TLR9
Species	Human
GeneID (Human)	54106
GeneID	54106
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 414-543 of human TLR9 (NP_059138.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR96
Protein Weight	116kDa
Shipping	Ice bag