

Anti-Tlx1 Polyclonal Antibody For Wb, Elisa - P31314

Item Number: CM0016043

Introduction

Rabbit polyclonal antibody against human TLX1, validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for spleen development studies and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	TLX1 Rabbit pAb
Remarks/Alias	TCL3; HOX11; TLX1
Species	Human
Gene ID	3195
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human TLX1 (NP_001182446.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31314
Protein Weight	34kDa
Shipping	Ice bag