

Anti-Tm7Sf2 Rabbit Polyclonal Antibody For Wb, Elisa - 076062

Item Number: CM0022280

Introduction

Anti-TM7SF2 rabbit polyclonal antibody suitable for WB and ELISA applications. Reacts with human and mouse TM7SF2, a delta(14)-sterol reductase involved in cholesterol biosynthesis. Ideal for lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TM7SF2 Rabbit pAb
Remarks/Alias	ANG1; C14SR; NET47; DHCR14A; TM7SF2
Species	Human
GeneID (Human)	7108
GeneID	7108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-250 of human TM7SF2 (NP_003264.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O76062
Protein Weight	46kDa
Shipping	Ice bag