

Anti-Tmc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8TDi8

Item Number: CM0027354



Introduction

High-quality rabbit polyclonal antibody against TMC1, a mechanotransduction channel subunit. Validated for WB, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for auditory system research.

[Learn More](#)

Parameter	Value
Product Name	TMC1 Rabbit pAb
Remarks/Alias	DFNB7; DFNA36; DFNB11; TMC1
Species	Human
GeneID (human)	117531
GeneID	117531
Immunogen	Synthetic peptide, a synthetic peptide corresponding to a sequence within amino acids 600-700 of human TMC1 (NP_619636.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDI8
Protein Weight	88kDa
Shipping	Ice bag