

Anti-Tmc5 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Uxy8

Item Number: CM0021914

Introduction

Rabbit polyclonal antibody targeting human TMC5 (Transmembrane channel-like protein 5). Validated for WB and ELISA; detects a probable ion channel component. Suitable for membrane protein research.

[Learn More](#)

Parameter	Value
Product Name	TMC5 Rabbit pAb
Remarks/Alias	TMC5
Species	Human
GeneID (human)	79838
GeneID	79838
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 676-760 of human TMC5 (NP_079056.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6UXY8
Protein Weight	115kDa
Shipping	Ice bag