

Anti-Tmed3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y3Q3

Item Number: CM0024658

Introduction

Rabbit polyclonal antibody targeting human TMED3, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studies of vesicular trafficking, Golgi dynamics, and early secretory pathway function.

[Learn More](#)

Parameter	Value
Product Name	TMED3 Rabbit pAb
Remarks/Alias	p26; P24B; p24g4; C15orf22; p24gamma4; TMED3
Species	Human
GeneID	23423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-180 of human TMED3 (NP_031390.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3Q3
Protein Weight	25kDa
Shipping	Ice bag