


Anti-Tmed9 Polyclonal Antibody For Wb And Elisa - Q9Bvk6

Item Number: CM0026699



Introduction

TMED9 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat TMED9. Targeting the p24 family protein involved in vesicular trafficking and Golgi retention, suitable for secretory pathway research.

[Learn More](#)

Parameter	Value
Product Name	TMED9 Rabbit pAb
Remarks/Alias	p25; GMP25; p24a2; HSGP25L2G; p24alpha2; TMED9
Species	Human
Gene ID	54732 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human TMED9 (NP_059980.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BVK6
Protein Weight	27kDa
Shipping	Ice bag