

Anti-Tmeff2 Monoclonal Antibody For Wb, Elisa - Q9Uik5

Item Number: CM0012127

Introduction

High-specificity TMEFF2 Rabbit Monoclonal Antibody validated for WB and ELISA. Reacts with Human and Rat TMEFF2 (tomoregulin-2). Ideal for neuroscience and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	TMEFF2 Rabbit mAb
Remarks/Alias	TR; HPP1; TPEF; TR-2; TENB2; CT120.2; TMEFF2
Species	Human
Gene ID (Human)	23671
Gene ID	23671
Immunogen	Synthetic peptide corresponding to aa 275-374 of human TMEFF2 (Q9UIK5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UIK5
Protein Weight	41kDa
Shipping	Ice bag