

Anti-Tmem100 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nv29

Item Number: CM0020424

Introduction

High-quality rabbit polyclonal antibody against human TMEM100, suitable for WB and ELISA. Cross-reacts with mouse. Validated for research use in vascular biology and pain signaling studies. 14 kDa transmembrane protein.

[Learn More](#)

Parameter	Value
Product Name	TMEM100 Rabbit pAb
Remarks/Alias	TMEM100
Species	Human
GenelD (Human)	55273
GenelD	55273
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human TMEM100 (NP_001093110.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NV29
Protein Weight	14kDa
Shipping	Ice bag