

Anti-Tmem119 Monoclonal Antibody For Wb, If-P, Ihc-P, Mihc, Elisa - Q4V9L6

Item Number: CM0006833

Introduction

Rabbit monoclonal antibody against human TMEM119 (OBIF), validated for WB, IF-P, IHC-P, mIHC, and ELISA. Cross-reacts with mouse and rat. Suitable for bone formation, microglia, and spermatogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	TMEM119 Rabbit mAb
Remarks/Alias	OBIF
Species	Human
GenelD (human)	338773
GenelD	338773
Immunogen	Recombinant protein of human TMEM119
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q4V9L6
Protein Weight	29kDa
Shipping	Ice bag