

Anti-Tmem119 Polyclonal Antibody For Wb, Elisa - Q4V9L6

Item Number: CM0022928

Introduction

Rabbit polyclonal antibody against human TMEM119 (OBIF), suitable for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. TMEM119 functions in osteoblast differentiation, bone mineralization, and is a microglia marker. Predicted molecular weight ~29 kDa.

[Learn More](#)

Parameter	Value
Product Name	TMEM119 Rabbit pAb
Remarks/Alias	OBIF; TMEM119
Species	Human
Gene ID (Human)	338773
Gene ID	338773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-270 of human TMEM119 (NP_859075.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q4V9L6
Protein Weight	29kDa
Shipping	Ice bag