

Anti-Tmem132D Rabbit Polyclonal Antibody For Wb, Elisa - Q14C87

Item Number: CM0015866

Introduction

TMEM132D Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat TMEM132D, a neuronal morphology regulator. Ideal for neuroscience research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	TMEM132D Rabbit pAb
Remarks/Alias	TMEM132D; MOLT; PPP1R153
Species	Human
GeneID	121256
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 730-915 of human TMEM132D (NP_597705.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14C87
Protein Weight	122kDa
Shipping	Ice bag