


Anti-Tmem147 Polyclonal Antibody For Wb, Elisa - Q9Bvk8

Item Number: CM0019082



Introduction

Rabbit polyclonal antibody targeting human TMEM147, validated for Western blot and ELISA. Ideal for studying multi-pass translocon complex function and endoplasmic reticulum membrane protein insertion. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	TMEM147 Rabbit pAb
Remarks/Alias	NEDFLPH; NIFIE14; TMEM147
Species	Human
Gene ID (Human)	10430
Gene ID	10430
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TMEM147 (NP_116024.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BVK8
Protein Weight	25kDa
Shipping	Ice bag