

Anti-Tmem176B Rabbit Polyclonal Antibody For Wb, Elisa - Q3Ybm2

Item Number: CM0023799

Introduction

Rabbit polyclonal antibody targeting TMEM176B (LR8), validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for dendritic cell maturation and cerebellar development studies. – Q3YBM2

[Learn More](#)

Parameter	Value
Product Name	TMEM176B Rabbit pAb
Remarks/Alias	TMEM176B; LR8; MS4B2; transmembrane protein 176B
Species	Human
Gene ID	28959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human TMEM176B (NP_001094781.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q3YBM2
Protein Weight	25kDa/29kDa
Shipping	Ice bag