

Anti-Tmem176B Rabbit Polyclonal Antibody For Wb, Elisa - Q3Ybm2

Item Number: CM0018922

Introduction

Rabbit polyclonal antibody against human TMEM176B (LR8, MS4B2). Validated for WB and ELISA; reacts with human and mouse. Suitable for dendritic cell maturation and cerebellar granule cell development research.

[Learn More](#)

Parameter	Value
Product Name	TMEM176B Rabbit pAb
Remarks/Alias	LR8; MS4B2; TMEM176B
Species	Human
Gene ID	28959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 148-208 of human TMEM176B (NP_054739.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q3YBM2
Protein Weight	29kDa
Shipping	Ice bag