

Anti-Tmem18 Polyclonal Antibody For Wb, Elisa - Q96B42

Item Number: CM0023418

Introduction

TMEM18 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and rat samples. Suitable for studying transcription repression, cell migration, and DNA-binding functions.

[Learn More](#)

Parameter	Value
Product Name	TMEM18 Rabbit pAb
Remarks/Alias	IncND; TMEM18
Species	Human
GeneID (human)	129787
GeneID	129787
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 40-100 of human TMEM18 (NP_690047.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96B42
Protein Weight	16kDa
Shipping	Ice bag