

Anti-Tmem192 Monoclonal Antibody For Wb, Elisa - Q8Iy95

Item Number: CM0002736

Introduction

Recombinant rabbit monoclonal antibody targeting human TMEM192, validated for Western blot (WB) and ELISA. Reacts with human and rat samples. Ideal for studying lysosomal and late endosomal membrane proteins.

[Learn More](#)

Parameter	Value
Product Name	TMEM192 Rabbit mAb
Species	Human
GeneID	201931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-271 of human TMEM192 (NP_001093859.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IY95
Protein Weight	31kDa
Shipping	Ice bag