


Anti-Tmem2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uhn6

Item Number: CM0020491



Introduction

Rabbit polyclonal antibody targeting human TMEM2 (CEMIP2), a cell surface hyaluronidase. Validated for WB and ELISA, cross-reactive with mouse and rat. Supports hyaluronan metabolism and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	TMEM2 Rabbit pAb
Remarks/Alias	TMEM2
Species	Human
Gene ID (Human)	23670
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1150-1250 of human TMEM2 (NP_037522.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UHN6
Protein Weight	154kDa
Shipping	Ice bag