

Anti-Tmem25 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q86Yd3

Item Number: CM0016056

Introduction

High-quality rabbit polyclonal antibody targeting human TMEM25. Validated for WB, IF/ICC, and ELISA. Detects TMEM25 in human, mouse, and rat samples. Ideal for studies on neuronal excitability and NMDA receptor degradation.

[Learn More](#)

Parameter	Value
Product Name	TMEM25 Rabbit pAb
Remarks/Alias	TMEM25
Species	Human
GeneID (Human)	84866
GeneID	84866
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 223-322 of human TMEM25 (NP_001137506.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YD3
Protein Weight	39kDa
Shipping	Ice bag