

Anti-Tmem43 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Btv4

Item Number: CM0003887

Introduction

Rabbit monoclonal antibody against TMEM43 (LUMA), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for nuclear envelope, innate immunity, and hearing research.

[Learn More](#)

Parameter	Value
Product Name	TMEM43 Rabbit mAb
Remarks/Alias	LUMA; ARVC5; ARVD5; AUNA3; EDMD7; AUNA2; TMEM43
Species	Human
Gene ID (Human)	79188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-135 of human TMEM43 (NP_077310.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTV4
Protein Weight	45kDa
Shipping	Ice bag