

Anti-Tmem61 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N0U2

Item Number: CM0021325

Introduction

Rabbit polyclonal antibody targeting human TMEM61, validated for Western blot (WB) and ELISA. Detects a 22 kDa transmembrane protein, ideal for membrane protein research and assay development.

[Learn More](#)

Parameter	Specification
Product Name	TMEM61 Rabbit pAb
Remarks/Alias	TMEM61
Species	Human
GeneID (Human)	199964
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 91-210 of human TMEM61 (NP_872338.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N0U2
Protein Weight	22kDa
Shipping	Ice bag