

Anti-Tmf1 Rabbit Polyclonal Antibody For Wb, Elisa - P82094

Item Number: CM0024436

Introduction

High-quality rabbit polyclonal antibody targeting human TMF1 (ARA160). Validated for WB and ELISA, cross-reacts with mouse. Ideal for transcription regulation and retrograde transport research.

[Learn More](#)

Parameter	Value
Product Name	TMF1 Rabbit pAb
Remarks/Alias	TMF; ARA160; TMF1
Species	Human
GeneID	7110
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TMF1 (NP_009045.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P82094
Protein Weight	123kDa
Shipping	Ice bag