

Anti-Tmlhe Rabbit Polyclonal Antibody For If/Icc And Elisa - Q9Nvh6

Item Number: CM0025352

Introduction

Rabbit polyclonal antibody against human TMLHE, a mitochondrial enzyme in carnitine biosynthesis. Validated for IF/ICC and ELISA. Targets trimethyllysine dioxygenase; relevant for metabolism and autism research.

[Learn More](#)

Parameter	Value
Product Name	TMLHE Rabbit pAb
Remarks/Alias	TMLD; TMLH; BBOX2; AUTSX6; TMLHED; XAP130; TMLHE
Species	Human
Gene ID (Human)	55217
Gene ID	55217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 167-376 of human TMLHE (NP_001171726.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NVH6
Protein Weight	50kDa
Shipping	Ice bag