

## Anti-Tmod1 Polyclonal Antibody For Wb And Elisa - P28289

Item Number: CM0027454



### Introduction

TMOD1 Rabbit pAb (P28289) is a polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat TMOD1. Tropomodulin-1 caps actin filament pointed ends. Ideal for cytoskeleton and actin dynamics research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | TMOD1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | TMOD; ETMOD; D9S57E; TMOD1                                                                                       |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 7111                                                                                                             |
| GeneID           | 7111                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-359 of human TMOD1 (NP_003266.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P28289                                                                                                           |
| Protein Weight   | 41kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |