

Anti-Tmod4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nzq9

Item Number: CM0027544

Introduction

Rabbit polyclonal antibody targeting human TMOD4, validated for WB and ELISA. Cross-reacts with mouse and rat. TMOD4 is a skeletal muscle tropomodulin that caps actin filaments, regulating cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	TMOD4 Rabbit pAb
Remarks/Alias	SK-TMOD; TMOD4
Species	Human
GeneID (human)	29765
GeneID	29765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-345 of human TMOD4 (NP_037485.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NZQ9
Protein Weight	39kDa
Shipping	Ice bag