


Anti-Tmpo Rabbit Polyclonal Antibody - P42166 / P42167

Item Number: CM0020610



Introduction

High-quality rabbit polyclonal antibody targeting human TMPO (LAP2). Suitable for WB, IHC-P, IF/ICC, and ELISA. Exhibits cross-reactivity with mouse and rat. Detects multiple isoforms (38-75 kDa). Ideal for nuclear architecture and immunology studies.

[Learn More](#)

Parameter	Specification
Product Name	TMPO Rabbit pAb
Remarks/Alias	TP; LAP2; CMD1T; LEMD4; PRO0868; TMPO
Species	Human
Gene ID (Human)	7112
Gene ID	7112
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-410 of human TMPO (NP_001027454.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42166/P42167
Protein Weight	38kDa/50kDa/75kDa
Shipping	Ice bag