

Anti-Tmprss11A Polyclonal Antibody For Wb And Elisa - Q6Zmr5

Item Number: CM0025350



Introduction

Rabbit polyclonal antibody against TMPRSS11A, a serine protease involved in cellular senescence. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TMPRSS11A Rabbit pAb
Remarks/Alias	HESP; ECRG1; HATL1; TMPRSS11A
Species	Human
Gene ID	339967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-421 of human TMPRSS11A (NP_872412.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZMR5
Protein Weight	47kDa
Shipping	Ice bag