

Anti-Tmprss11A Polyclonal Antibody For Wb, Elisa - Q6Zmr5

Item Number: CM0013843



Introduction

Rabbit polyclonal antibody against human TMPRSS11A, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Useful for senescence and esophageal cancer research. Predicted MW 47 kDa.

[Learn More](#)

Parameter	Value
Product Name	TMPRSS11A Rabbit pAb
Remarks/Alias	HESP, ECRG1, HATL1, TMPRSS11A
Species	Human
GeneID (Human)	339967
GeneID	339967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-421 of human TMPRSS11A (NP_872412.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZMR5
Protein Weight	47kDa
Shipping	Ice bag