

Anti-Tmprss2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15393

Item Number: CM0020136



Introduction

TMPRSS2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human TMPRSS2 serine protease involved in viral entry and prostate cancer. Cross-reacts with mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	TMPRSS2 Rabbit pAb
Remarks/Alias	PRSS10; TMPRSS2
Species	Human
Gene ID (Human)	7113
Gene ID	7113
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TMPRSS2 (NP_005647.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15393
Protein Weight	54kDa
Shipping	Ice bag