

Anti-Tmprss2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O15393

Item Number: CM0008949



Introduction

High-quality rabbit monoclonal antibody against human TMPRSS2 (PRSS10) for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Supports research on host-virus interactions, cancer, and protease biology.

[Learn More](#)

Parameter	Value
Product Name	TMPRSS2 Rabbit mAb
Remarks/Alias	PRSS10; TMPRSS2
Species	Human
GeneID (Human)	7113
GeneID	7113
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TMPRSS2 (O15393).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15393
Protein Weight	54kDa
Shipping	Ice bag