

Anti-Tmprss6 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q8IU80

Item Number: CM0013743



Introduction

Anti-TMPRSS6 Rabbit Polyclonal Antibody validated for WB, IF-P, ELISA in human, mouse, rat. Targets matriptase-2 for iron homeostasis research. Recombinant protein immunogen, 90 kDa.

[Learn More](#)

Parameter	Specification
Product Name	TMPRSS6 Rabbit pAb
Remarks/Alias	MT2; IRIDA; TMPRSS6
Species	Human
GeneID (Human)	164656
GeneID	164656
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 77-220 of human TMPRSS6, NP_705837.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IU80
Protein Weight	90kDa
Shipping	Ice bag