

Anti-Tmtc3 Polyclonal Antibody For Wb And Elisa - Q6Zxv5

Item Number: CM0017781

Introduction

Rabbit polyclonal antibody against human TMTC3 (O-mannosyl-transferase), suited for WB and ELISA. Immunogen corresponds to amino acids 615-914. Useful for investigating ER stress response and cadherin O-mannosylation.

[Learn More](#)

Parameter	Value
Product Name	TMTC3 Rabbit pAb
Remarks/Alias	LIS8; SMILE; TMTC3
Species	Human
GeneID (Human)	160418
GeneID	160418
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 615-914 of human TMTC3 (NP_861448.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6ZXV5
Protein Weight	104kDa
Shipping	Ice bag