


Anti-Tmx1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H3N1

Item Number: CM0018428



Introduction

Anti-TMX1 polyclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, rat TMX1 protein (32kDa). Ideal for redox biology and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	TMX1 Rabbit pAb
Remarks/Alias	TMX; TXNDC; PDIA11; TXNDC1; TMX1
Species	Human
GeneID	81542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-280 of human TMX1 (NP_110382.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3N1
Protein Weight	32kDa
Shipping	Ice bag