


Anti-Tmx1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H3N1

Item Number: CM0019614



Introduction

Rabbit polyclonal antibody targeting human TMX1, a thioredoxin-related transmembrane protein involved in protein disulfide isomerase activity, ER stress, and triglyceride metabolism. Validated for Western blot (WB) and ELISA.

[Learn More](#)

Parameter	Value
Product Name	TMX1 Rabbit pAb
Remarks/Alias	TMX; TXNDC; PDIA11; TXNDC1; TMX1
Species	Human
GeneID (Human)	81542
GeneID	81542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-280 of human TMX1 (NP_110382.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H3N1
Protein Weight	32kDa
Shipping	Ice bag