


Anti-Tnf-A Monoclonal Antibody For Ihc-P, Elisa - P01375

Item Number: CM0007398



Introduction

Rabbit monoclonal antibody specific to human TNF- α (cachectin, TNFSF2) with validated IHC-P and ELISA applications. Immunogen: recombinant fusion protein covering aa 77-233. Ideal for studying inflammation, cancer, and autoimmune disorders.

[Learn More](#)

Parameter	Value
Product Name	TNF- α Rabbit mAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha; TNF- α
Species	Human
Gene ID (Human)	7124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-233 of human TNF- α (NP_000585.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01375
Protein Weight	26kDa
Shipping	Ice bag