

Anti-Tnf-A Rabbit Monoclonal Antibody For Wb, Elisa - P06804

Item Number: CM0006578



Introduction

Mouse TNF- α rabbit monoclonal antibody validated for WB and ELISA, recognizes both mouse and rat TNF- α . Ideal for research in inflammation, immune response, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	TNF- α Rabbit mAb
Remarks/Alias	DIF; Tnfa; TNF-a; TNFSF2; Tnlg1f; Tnfsf1a; TNFalpha; TNF-alpha; TNF- α
Species	Mouse
GeneID (Human)	7124
GeneID	21926
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-235 of mouse TNF- α (NP_038721.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P06804
Protein Weight	26kDa
Shipping	Ice bag